





Abstract

- The US Food and Drug Administration (FDA) recently released a new draft guidance enabling generic manufacturers of peptide drugs to file an Abbreviated New Drug Application (ANDA) for synthetic peptide drug products.
- Processes for manufacturing the generic and reference drug (RLD) are not equivalent, leading to manufacturing related impurities.
- Manufacturers are required to prove that the synthetic peptide product does not contain impurities with an increased risk of immunogenicity that could result in the development of anti-drug antibodies.
- We use both in silico analysis and in vitro validation assays to perform immunogenicity risk assessment of peptide generics. This process is referred to as the **PANDA assay** which can be used to support generic peptide drug equivalency in an ANDA application.


PANDA Overview



Impact of Manufacturing Impurities on T cell Epitopes




HLA Class II Binding Assay



HLA-Class II Binding Assay: Peptides that are predicted to bind HLA are synthesized and assayed over a range of 7 concentrations in our HLA-binding assay. In brief, peptides are incubated overnight with soluble HLA and a biotin labeled competitor of moderate affinity. On day 2, the reaction is halted and the mixture is transferred to a plate coated with a pan anti-HLA antibody. On day 3, plates are developed by the addition of streptavidin-Europium and fluorescence is measured.


Epitope Prediction by EpiMatrix



In Vitro Immunogenicity Protocol (IVIP)

- The ability of the test article (new Generic and impurities) and the RLD to stimulate a *de novo* T-cell response is compared to several controls including HSA (protein neg control), KLH (protein positive control) and a CEFT (protein pool positive control).
- 14 days post exposure, cells are harvested and plated into precoated IFN γ ELISpot plates. Cells are restimulated and incubated overnight. On day 15, ELISpot plates are developed and sent to Zellnet Consulting Inc. for blind, independent analysis.

Timeline of Naive Immune Response



Summary of IFN γ Fluorospot responses across donors – RLD A

DONOR	1	2	3	4	5	6	7	8	9	10	TOTAL # of Positive Responses*
DRB1 Allele	04:03 15:01	01:01 07:01	07:01 07:01	03:01 16:01†	03:01 07:01	01:01 14:01†	03:01 13:03†	03:01 09:01†	01:02 07:01	03:01 15:01	
EpiMatrix Hits: Allele	1 0	1 1	1 1	1 0	1 1	1 NA	1 NA	1 0	1 1	1 0	
RLD - A	1	-	-	-	-	-	-	-	-	-	0/10
	2	+	-	+	-	-	-	-	-	-	2/10
Test Article - A	11	-	-	-	-	-	-	+	-	+	2/10
	12	-	+	-	-	+	-	-	-	-	2/10
	13	-	+	-	-	-	-	-	-	-	1/10

* A response is considered positive with >50 SFC/10⁶ cells and a stimulation index ≥ 2
† Epitope predictions (EpiMatrix Hits) for HLA DRB1*16:01 were modeled on supertype relative DRB1*1101; predictions for HLA DBR1*14:01 and *13:03 are not available; predictions for HLA DRB1*0901 are available through EpiVax internal models not normally included in PreDeFT analysis.

References

- Steere AC, Kitz W, Drouin EE, et al. Antibiotic-refractory Lyme arthritis is associated with HLA-DR molecules that bind a Borrelia burgdorferi peptide. *J Exp Med* 2006; 203: 961–971.
- Wullner D, Zhou L, Bramhall E, et al. Considerations for optimization and validation of an in vitro PBMC derived T cell assay for immunogenicity prediction of biotherapeutics. *Clin Immunol* 2010; 137: 5–14.
- *Lund et al. Definition of Supertypes for HLA Molecules Using Clustering of Specificity Matrices. *Immunogenetics*. 2004; 55(12):797–810.
- **Southwood et al. Several Common HLA-DR Types Share Largely Overlapping Peptide Binding Repertoires. *J Immunol*. 1998; 160(7):3363–73.



Modeling Unnatural Amino Acids

- When EpiMatrix cannot model HLA binding for modifications found in the impurity, we use a sensitivity analysis to find the best proxy.
- The modified residue is replaced with a neutral placeholder "X". We then replace "X" with each of the natural 20 amino acids.
- The goal is to determine if any residue at these positions can lead to a significant increase or decrease in predicted HLA binding potential.
- We also compare the properties of the chemically modified residues with the naturally occurring amino acids and pick the "best-matched" residue as a proxy

